





































# Bioconnect SpO<sub>2</sub> Sensors







































FEATURES	BENEFITS
Ergonomic design	Comfortable fit. Proper finger placement for accurate readings
Spring is integral part of clip	Can not pop out and maintains proper tension
Stainless steel axle	Will not break. Adds to overall strength of product
Vari-torque design of strain relief	Reduces cable breakage
Cable anchored at both ends with Kevlar	Stronger, increasing life span of cable
Shielded cable	Truer, cleaner signal
PVC free polyurethane jacket	Environmentally sensitive and patient friendly
Clearly labeled packaging	Ease in identification in hospital
Instructions in every package	Ensures proper usage
Cable labels	Prevents unwanted usage
One year warranty	Confidence in product
Priced competitively	Substantial savings over OEMs
Compatible with most OEMs	Wider applications
Reusable	More cost efficient
Repairable	Long term financial savings

	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ext ension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Agilent Technologies Monitors</b>								
	<b>SEE PHILIPS / HEWLETT PACKARD</b>							
<b>Alaris (VAC)</b>								
Vital Check 4400 Series	1002N	4'/1.2m	1002E	10'/3m	5451; Extension: 5459			
	1002N-10	10'/3m	N/A	-	511-10L			
<b>Allegiance Healthcare</b>								
Allegiance	<b>SEE MASIMO</b>							
<b>Alliance Instruments</b>								
Integriti Model	<b>SEE NOVAMETRIX (1016)</b>							
<b>AMI</b>								
Nellcor/Invivo Based (Apnea systems may use standard Nellcor adapters or interface cable)	1011	4'/1.2m	1036A	10'/3m	Sensor: DS-100A			
	1036	10'/3m	N/A	-				
Criticare based System 7310 / 7500	1002N	4'/1.2m	1002A	10'/3m	5451; Extension: 5459			
	1002	10'/3m			511-10L			
<b>Artema</b>								
ATHENA 9060 (OHMEDA Technology)	1021	10'/3m	N/A	-	Sensor: SL033057		Uses OxyTip® interconnects Direct Connect Available Only	
















	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Agilent Technologies Monitors</b>								
	<b>SEE PHILIPS / HEWLETT PACKARD</b>							
<b>Alaris (VAC)</b>								
Vital Check 4400 Series	1002N	4'/1.2m	1002E	10'/3m	5451; Extension: 5459			
	1002N-10	10'/3m	N/A	-	511-10L			
<b>Allegiance Healthcare</b>								
Allegiance	<b>SEE MASIMO</b>							
<b>Alliance Instruments</b>								
Integriti Model	<b>SEE NOVAMETRIX (1016)</b>							
<b>AMI</b>								
Nelcor/Invivo Based (Apnea systems may use standard Nelcor adapters or interface cable)	1011	4'/1.2m	1036A	8'/2.4m	Sensor: DS-100A			
	1036	10'/3m	N/A	-				
Criticare based System 7310 / 7500	1002N	4'/1.2m	1002A	8'/2.4m	5451; Extension: 5459			
	1002	10'/3m			511-10L			
<b>Armstrong</b>								
AD-900, 1000, 2000, 3000	1001	4'/3m	1011E	10'/3m	BCI			

	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Artema</b>								
ATHENA 9050, 9060 (OHMEDA Technology)	1021	10'/3m	N/A	-	Sensor: SL033057		Uses OxyTip® interconnects Direct Connect Available Only	
<b>Artema Continued</b>								
Gamme 8000, Diascope 2, S&W (new), ATHENA (new) (NELLCOR technology)	1011	4'/1.2m	1021A	8'/2.4m	Sensor: DS-100A			
	1021B	10'/3m	N/A	-	N/A			
S&W (OLD), MM200	1050	10'/3m	N/A	-	Sensor: 50-70259-00			
<b>Atom Medical</b>								
<b>SEE MASIMO (1049)</b>								
<b>Baxter Medical</b>								
ASAT (Older Simed technology)	1031	10'/3m	N/A	-	Sensor: FP-100 Edwards C.C. Division Simed : 600015			
Explorer (Nelcor) / Monarch	1011	4'/1.2m	1031A	8'/2.4m				












	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>BCI</b>								
3100 series, 3300-3304, 9000 series, 9100, Advisor 9200, Autocorr Plus 3404, Mini-torr, Mini-torr Plus 6004	1001	4'/1.2m	1001E	10'/3m	Sensor: 3044 Sensor: 3444			
SERENA 6100 Series	1001	4'/1.2m	1001A	8'/2.4m	Sensor: 3044 Sensor: 3444			
<b>Bionet</b>								
BM3 (Nellcor Technology)	1011	4'/1.2m	1011E	10'/3m	Sensor: DS-100A Extension: EC-8 (optional)			
	1011N	10'/3m	N/A					
<b>Biosys</b>								
OLD GENERATION BPM-300	1069	10'/3m	N/A	-				
NELLCOR TECH. 700, Guardian BPM-200	1069N	4'/1.2m	1011E	10'/3m	Sensor: BSR-250A Extension: EC-8 (optional)			
<b>Bird</b>								
4400 (CRITICARE TECHNOLOGY)	1002N	4'/1.2m	1002A	8'/2.4m	Sensor Direct: 511-10L Sensor (4'): 511-SD			
	1002	10'/3m	N/A		Adapter: 5459A			
<b>Bruker</b>								
<b>ODAM (BCI Tech.)</b> 40,493,444	1001	4'/1.2m	1075C	8'/2.4m	Sensor : 3044 Sensor : 3444			

	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Burdick</b>								
<b>OXY100, 200, PRO2 (BCI Tech.)</b>	1001	4'/1.2m	1001A	8'/2.4m	Sensor : 3044 Sensor : 3444			
<b>Charter</b>								
	<b>SEE KONTRON</b>							
<b>Colin</b>								
BP88/306/508/608/808,885 PressmateBX10ne (Nellcor Tech.) (NOTE:BX308/408/BX5 SENSORS UNAVAILABLE) NEWER OXI-MAX Sensors for Prodigy Systems Unavailable - Call for Repair info	1011	4'/1.2m	1011E	10'/3m	Sensor: DS-100A Extension: EC-8 (optional)			
BP 88S (Nellcor Tech.) BP308,408,BXS	1011	4'/1.2m	1011A	8'/2.4m	Sensor: DS-100A Adapter: SCP-10 / MC-10			
	1011P	10'/3m	N/A	-				
BP308,408, BX5, BX10ma (Masimo Tech.)	N/A - repair only							
<b>Corometrics</b>								
118/188/120/128/129/510/511/ 556 (Nellcor Tech.) Newer GE Corometrics may require 1033 connector (SEE GE)	1011	4'/1.2m	1040A	8'/2.4m	Sensor: DS-100A Adapter: 4033CAX			
	1040	10'/3m	N/A	-				
511 (Ohmeda Tech.)	1017	10'/3m	N/A	-	6051-0000-035			
























 Bioconnect Division of RF Industries	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
Corpus								
08/16 (BCI Tech.) (Light Blue Conn.)	1001	4'/1.2m	see note		Unit requires use of OEM adapter cable. Light Blue Color (CALL FOR INFO)			
08/16 08/16S (Masimo Tech.) (Red Conn.)	N/A - repair only							
Criticare (CSI)								
501/503/504/507/508/581 1100 POET/ 1100PONI/ Scholar (Spot, Plus, and DX models may NOT need adapter Use 1002E for Extension or 1002N-10)	1002N	4'/1.2m	1002A	8'/2.4m	Sensor: 511-SD Direct Connect: 511-10L			
	1002N-10	10'/3m	N/A	-				
	1002	10'/3m	N/A	-	511-10L			
Critikon(GE) Sensormedics Dinamap								
Oxyshuttle Old Generation	1043S	3'/.9m	1003A	10'/3m	Sensor: 9084 Adapter: 9025			
OxyTrak, 8280, 1846SX (w/8802 adapter)	1003N	4'/1.2m	see note		Sensor: 8997 Adapter.: 8802 (REQUIRED) N/A			
DINAMAP + (OLD Generation) 8700 Series Oxyshuttle New Generation	1043S	3'/.9m	1027A	10'/3m	Sensor: 9084 or 766721 (Sensormedics) Adapter: 8731			
	1003P	10'/3m	N/A	-				




















	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Critikon(GE) Sensormedics Dinamap</b>	<b>Continued</b>							
DINAMAP + (NEW Generation) 9700 Series MPS, SELECT, PORTABLE & COMPACT Pro Care (Nelcor Tech.) Some OxiMax Pro Care require 3M 14-P adapter	1011	4'/1.2m	1011E	10'/3m	Sensor: DS-100A Ext.: EC-8 (Optional) NOTE: May require DOC-10 adapter			
PRO CARE 400 (Masimo Tech.)	N/A - repair only							
<b>Datascope</b>								
Accusat, Accutor 3/4/SAT/PLUS, Passport (Datascope Tech.)	1004N	4'/1.2m	1004A	8'/2.4m	0600-00-0026-01			
	1004	10'/3m	N/A	-	0600-00-0026-02			
Passport II/IILT, EL, XG Accutor Plus, Expert, (Nelcor tech.) (Trio/Spectrum w/opt. Nelcor)	1011	4'/1.2m	1011A	8'/2.4m	DS-100A			
	1011P	10'/3m	N/A	-	MC-10 or SCP-10			
Passport II/IILT, XG Accutor Plus, Trio, Spectrum (Masimo Tech.)	N/A - repair only							
<b>Datex- ENGSTROM (DATEX/OHMEDA)-GE</b>								










**SEE OHMEDA**

















	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Digicare</b>								
BCI Technology	1001	4'/1.2m	1051A	8'/2.4m	Sensor: 3044 Adapter: CP599-0195			
	1051	10'/3m	N/A	-				
<b>Dimi Italia</b>								
	<b>SEE NELLCOR</b>							
<b>Dinamap</b>								
	<b>SEE CRITIKON</b>							
<b>Dixtal</b>								
	<b>SEE NOVAMETRIX (1016)</b>							
<b>Draeger</b>								
Oxipac 2000, 2500 (BCI Tech.)	1001	4'/1.2m	1001E	10'/3m	Sensor: 3044 1001E Extension Optional			
Vitalert 1000, Cicero, Narcormed, Cato, Vitaray PM8014, 8060 (Nellcor Tech.)	1011	4'/1.2m	1070A	8'/2.4m	Sensor: DS-100A Adapter: 5720084			
	1070	10'/3m	N/A	-				
OxySat 2, NarkoMed (Nellcor Tech.)	1011	4'/1.2m	1011PA	13'/4m	Sensor: DS-100A Adapter: Pre-Amp			
Newer Generation Siemens DragerMedical Monitors, Infinity Series	<b>SEE SIEMENS</b>							

	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Elmed</b>								
<b>PSO1000</b>	<b>SEE B.C.I. (1001)</b> Sensor: 6006194							
<b>Entri Medical</b>								
ProPac 102/104/106 (Nellcor Tech.)	<b>1011</b>	4'/1.2m	1011PA	13'/4,	Sensor: DS-100A Adapter: Pre-Amp			
<b>Fukuda Denshi/ Dynascope</b>								
HG302 MODULE/ DS Series	<b>SEE NELLCOR 1011</b>							
<b>Gambro</b>								
EOS (Criticare Tech.) Some monitors (BP806) may use Nellcor tech.	<b>SEE CRITICARE</b> (Requires Direct Connect 1002 or 1002A adapter) Sensor 511-10L							
<b>Getemed</b>								
VitalGaurd Nellcor or Masimo Based - Customer model specific	<b>SEE MASIMO (1049) or NELLCOR (1011)</b> Call for help or information							
<b>G.E. Medical Systems</b>								
Marquette: 7000 Series/TRAM A.R./TRAM Modular x00SL Series Hellige: SMU610 SMK611SMV210(Ohmeda Tech.)	<b>SEE OHMEDA (1017)</b>							
Marquette: Handhelds (Old Tech.)	<b>SEE NONIN (1014)</b>							
Marquette: DASH/Eagle/Solar/ TRAM x50SL/ TramScope/ Responder 3000 (Nellcor tech.)	1011	4'/1.2M	1033A  3300 Oximax	8'/2.4m	Sensor: DS-100A Adapter: D9004GE 2006644-001  (call for OxiMax E9004 1033OM)			
	1033B	10'/3m	N/A	-				
















	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
Marquette: DASH 3000/4000 Solar and Tramscope Modules with (Masimo tech.)	N/A for repair only							
FOR COROMETRIC CRITIKON DINAMAP SEE RESPECTIVE OEMs	<b>SEE RESPECTIVE OEMs</b>							
<b>Goldway</b>								
(UT6000A-F use Nellcor sensors)	<b>SEE NELLCOR 1011</b>							
UT4000F BCI	1001	4'/1.2m	8400	8'/2.4m				
	8400	10'/3m	N/A	-				
<b>Healthdyne/ Resperonics</b>								
Model 900 Series	<b>SEE NONIN (1014)</b> (Note Some Resperonics may use Novamterix S.B sensors (1016)) or Handhelds may use 1016N DB-9 with bottom key							
<b>Heska (VetOX, Vet CO)</b>								
4400 Series	<b>SEE B.C.I. (1001)</b>							
<b>Hewlett Packard</b>								
	<b>SEE PHILIPS</b>							
<b>Invacare</b>								
	<b>SEE B.C.I. (1001)</b>							
<b>Infinium</b>								
	<b>SEE NELLCOR 1011</b>							


















	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Invivo Research MDE</b>								
4500 Scout, Scout+/-1/2/3/4, Scout ES, 4500 Plus, Millenia, Omega 1445, Omni-Trak, TVS, Sat Scout	1007	10'3m	N/A	-				
Same as above (w / Nellcor tech.)	1011	4'/1.2m	1007A	8'/2.4m				
Millenia, 4500 (Masimo tech.)	N/A for repair only							
MDE, InvivoMDE, Escort (Masimo)	N/A for repair only							
Medical Device Equipment (M.D.E.) Escort Series Nellcor Based	<b>SEE M.D.E.</b>							
<b>Ivy Biomedical</b>								
405A, 405D, 405T, 405P, 705 (405A & 405D Require 1011PA Adapter)	1011	4'/1.2m	1011PA	13'/4m				
SatGaurd IVY 2000, 405TM, 405CM	Masimo - N/A for repair only							
<b>Kontron / Charter</b>								
Micro/Mini/Super/Kolormon, 7840, 7845, 7137(Nellcor Modules)	1011	4'/1.2m	1052A	8'/2.4m				
	1052	10'/3m	N/A	-				
Charter, VitalSat and (Masimo modules)	Masimo - N/A for repair only							




















	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Lynn Medical</b>								
	<b>SEE B.C.I. (1001)</b>							
<b>Magnetic Resonance Equip. Corp. (MRE)</b>								
MRI 9500 (ADAPTER ONLY)	N/a	-	1065A	8'2.4m				
<b>Masimo</b>								
All models	N/A - repair only							
<b>Medical Data Equipment (M.D.E.)</b>								
Escort I, Escort 100, OPT11A, Escort 300A	1011	4'/1.2m	1011PA	13'/4m	Sensor: DS-100A Adapter: C-13-200			
Escort II/OPT30/Prism/Type 20300/20401-03 (SEE InvivoMDE for more info)	1011	4'/1.2m	1011E		Sensor: DS-100A Extension: EC-8 (Optional)			
	1011N	10'/3m	N/A	-				
<b>Medical Equip. Korea (M.E.K.)</b>								
MP-100, 110, 400, 500, 600, 1000	7600	10'/3m	N/A					









	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Medical Research Lab (MRL)</b>								
ALL	1047	4'/1.2m	N/A	-				
<b>Medlab</b>								
MOST MONITORS (BCI Tech.)	1001	4'/1.2m	1001A	8'/2.4m				
	1001R	10'/3m	N/A	-				
MONITORS WITH DUAL CONN. KEYS (BCI tech.)	1001	4'/1.2m	1074A	8'/2.4m				
	1074	10'/3m	N/A	-				
<b>Mennen</b>								
HORIZON XL (OHMEDA Tech.)	<b>SEE OHMEDA (1017)</b>							
Horizon XL (Nellcor), Mercury 1200, MR 1300, Envoy	1011	4'/1.2m	1053A	8'/2.4m	Sensor: DS-100A Adapter: 261-877-070			
	1053	10'/3m	N/A	-				
Horizon 1100/2000XL Mercury	1011	4'/1.2m	1011PA	13'/4m	Sensor: DS-100A Adapter: Pre-Amp 1011PA (Required)			
























		BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
Mindray									
Criticare monitors MEC-509,509+,503,507,	1002N	4'/1.2m	1002A	8'/2.4m	Adapter 512A-30-06074				
	1002	10'/3m	N/A	-					
PM-6201,7000,8000,9000	1011	4'/1.2m	8300	8'/2.4m	Sensor DS-100A Adapter 0100-20-42594				
	8300	10'/3m	N/A	-					
Moriya									
SD Series (SD01-BC)		SEE B.C.I. (1001)							
MSA / Catalyst Research / Miniox									
Mini Ox IV	1042N	4'/1.2m	1042A	8'/2.4m					
	1042	10'/3m	N/A	-					
Mini Ox V	1042V	10'/3m	N/A	-					
Nellcor									
N-100 (4 Button) Call For info	1011	4'/1.2m	N/A	-	Sensor: DS-100A NOTE: PRE-AMP REQUIRED (NOT AVAILABLE)				


	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Nellcor</b>	<b>Continued</b>							
N-100 (5 Button) Hi Sat, N-200	1011	4'/1.2m	1011PA	13'/4m	Sensor: DS-100A Adapter: C-13-200			
N-180, NPB-40	1011	4'/1.2m	1011E	10'/3m	Sensor: DS-100A Extension: (optional) EC-8			
	1011N	10'/3m	N/A	-				
N-/190/290/295/350/390/395 N- 3000 N-6000 Symphony (NON OXIMAX)	1011	4'/1.2m	1011A	8'/2.4m	Sensor: DS-100A Adapter: MC-10 /SCP-10			
	1011P	10'/3m	N/A	-				
N65, NPB40 (upgraded to Oximax)	N/A for repair only	4'/1.2m						
N-550/595 OXIMAX Compatible	N/A for repair only	4'/1.2m	1011B	8'/2.4m	Call For Repair Info and Availability (MUST use OEM sensors) Adapter: DOC-10			
<b>Nihon Kohden (N.K.)</b>								
AL-900P Module, BSM- 1102/9510/4104A/4114A, OGS 2000 (Nellcor Technology)	1013N	4'/1.2m	1013A	8'/2.4m	TL-202 S			
	1013	10'/3m	N/A	-	TL-101S			

	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Nihon Kohden (N.K.) Continued</b>								
AL-900P Module, BSM-1101,1102/9510/4104A/4114A, OGS 2000 (Nellcor Technology)	1011	4'/1.2m	1011A 1301 (Oximax)	8'/2.4m	Sensor: DS-100A Adapter: MC-10 /SCP-10			
	1011P	10'/3m	N/A	-				
LifeScope (NEW) i/N/A/P/M/S BSS-9800/ BSM-4000/5000/9000/2300 OPV-1500	1300	10'/3m	N/A	-	SENSOR: TL-201T TL-220T, TL-631T3			
<b>Nonin</b>								
ALL (NOTE: XPOD Model 3102 adapter unrepairable)	1014	4'/1.2m	1014E	10'/3m	Sensor: 8000K2, AA,J Extension: (1014E optional)			
8604	1014	4'/1.2m	1014A	8'/2.4m	Sensor: 8000K2, AA,J Adapter: 8604D			
	1014L	10'/3m	N/A	-				
8600 (w/rear mount flat 8 pin conn.)	1014	4'/1.2m	1014B	8'/2.4m				
	1014C	10'/3m	N/A	-				
<b>Novamatrix</b>								
MODEL 500 (non-SuperBright) Also On Spacelabs 90419, 90501, 90651A-03	1015	10'/3m	1016E	10'/3m	Sensor: 8660 Spacelabs P/N: 369083-001 Extension Cable Optional		 Female	

















	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Novamatrix</b>	<b>Continued</b>							
505/515B- C/520/7000/Oxypleth/ MarSpO2 (Super Bright)	1016	10'/3m	N/A	-	Sensor: 8776-00 (OxySnap Conn. Unavailable) Spacelabs P/N: 015- 0130-00			
GE	1604	10'/3m	N/A	-				
Handheld Units 512/710/715 MODEL 510 (w/Oxysnap square conn. Directly to monitor NOT available)	1016N	4'/1.2m	N/A	-	Adapter for table top units available (0400-1600)			
<b>Odam (BCI Tech.)</b>								
BCI Technology.	<b>SEE BRUKER SEE ALSO SCHILLER</b>							
<b>Ohmeda (Datex Engstrom) GE</b>								
Ohmeda Model 2000, 405T (Masimo Technology) Call for Availability	Masimo - N/A for repair only							
3700, 3710, Biox 3740, Biox 3760, 4700, 3770, 3775, 3800, 3900 OxiCap, 5250RGM, Rascal II, TuffSat (Ohmeda Technology)	1017	10'/3m	N/A	-	Sensor: (OLD) 6051-0000-035 (New): Oxy-F4-H Mennen P/N: 800- 060-540			
AS-3, CG-25, CH-RS, CS-3, CapnoMan Ultima, CardioCap I, II, 5, M-Estpr, Oscar II, OSP- 200, ULTSVI/SVO/SIO/SOO, Satlite,	1005N	4'/1.2m	1005A	8'/2.4m	Sensor: (Old) SAS- F (New) Oxy-F-DB Direct: (Old) SAS- F4 (New) Oxy-F-N Adapter: Oxy-C3			
	1005	10'/3m						
OxyTip+ Sensors and Adapters	Circular OxyTip+ Adapters and Sensors unavailable. Direct connect or standard Sensors only. Call for info							














	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Ohmeda</b>	<b>Continued</b>							
3700, 3710, Biox 3740, Biox 3760, 4700, 3770, 3775, 3800, 3900 OxiCap, 5250RGM, Rascal II, TuffSat (Ohmeda Technology)	1700	4'/1.2m	0400-1700	8'/2.4m				
<b>Pace Tech</b>								
Model 520 SwitchCraft connector BCI Tech. VitalMax and Handy Pack use 1032N direct Use 1001E for Extension	1032N	4'/1.2m	1032A	8'/2.4m	Sensor: 4510 Adapter: 4530			
Mini-Pack 911-S, Newer Monitors, Redel Connectors. (BCI Tech.) Use B400-1001E for Extension on Handhelds	1032N	4'/1.2m	1032NA	8'/2.4m	Sensor: 4510 Adapter: 4530			
<b>Palco</b>								
ALL	1018	4'/1.2m	N/A	-	POX020-100			
	1018-10	10'/3m	N/A	-	POX020-105			
<b>Philips / Hewlett Packard (HP) / Agilent</b>								
(Nelcor Based Tech.) HP CMS,V24/V26 SpO2 module M1020A, Virdia Series 50XM/XMO/IP M1020A >=Rev C, HP CodeMaster XL/+100, 78352C, 78354C, M1205A, M1722A-B, M1723A-B, M2475B, Merlin M1020 Rev E or later; M1350C series 50XMO Fetal monitor	1011	4'/1.2m	1006A	8'/2.4m	Sensor: DS-100A/M1191T Adapter: M1900B			



















	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Philips / Hewlett Packard (HP) / Agilent</b>	<b>Continued</b>							
(Phillips Based Tech.) HP CMS,V24/V26 SpO2 module M1020A, Virdia Series 50XM/XMO/IP M1020A >= Rev C , HP CodeMaster XL/+100, 78352A/C, 78354A/C, 78833A- B, M1722A-B, M1723A-B, M2475B, Merlin M1020 Rev E or later; M1350 Series 50XMFetal monitor	1006V	4'/1.2m	0600	8'/2.4m	Sensor: M1191A Adapter: M1940A		  Female	
(Philips Based Technology) M2600A w/M2601A, MP30/30/40/50TeleMon A or B w/Transmitter A or B, C3:862474,862478, M3/M4: M3000A 3002 (Pot. A01 X2), Intellivue:, Module M1020B, HeartStream/Start XL/XLT, HeartStart 4000 HeartStart MRx; 8102 MP2; M8105 MP5; M2703A Avalon FM30 Fetal monitor; C3 (862474, 862478); Sure signs VM monitors (863063, 863064, 863065, 863066, 863068, 863077); Sure Signs VS3 monitor (863071, 863072, 863074)	1006V	4'/1.2m	0601	6'/2m	SensorM1191A  Extension M1941A Adapater: M1943A (1m) M1943AL (3m)		  Female	
	1006V-10	10'/3m	N/A	-	Sensor M1191A			





	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Quartz Medcial</b>								
All	Masimo - N/A for repair only							
<b>Resperonics</b>								
	<b>SEE HEALTHDYNE</b>							
<b>S&amp;W Vickers/ Artema</b>								
	<b>SEE ARTEMA</b>							
<b>Schiller</b>								
Argus LCM Plus (Masimo Tech.)	Masimo - N/A for repair only							
Oximeter OX-1, Argus Pro (Standard Masimo)	Masimo - N/A for repair only							
CM8 (BCI tech.)	1001	4'/1.2m	1075B	8'/2.4m				
Argus TM7 (Nellcor tech.)	1011	4'/1.2m	1075A	8'/2.4m				
ODAM BRUKER (BCI Tech.)	1001	4'/1.2m	1075C	8'/2.4m				
<b>Siemens</b>								
MircO2, SC6000, 7000, 8000, 9000, Sirecust 404-1, 700, 630, Sirem SpO2 Module (Nellcor Technology) NON-MultiMed Pods	1011	4'/1.2m	1024A	8'/2.4m	Sensor: DS-100A Adapter: 453384OEH50U			
	1024	10'/3m	N/A	-				
<b>Siemens / Dräger Medical</b>								
SC6000, SC7000, SC8000, SC9000 New Generation w/ MultiMed Pods	1011	4'/1.2m	7000	8'/2.4m	Sensor: DS-100A Adapter: 3368433			

	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Siemens / Dräger Medical</b>								
New Generation MultiMed Infinity Series SireCust or SC Series with Additional MultiMed Software,	MultiMed and Handheld unit have options for both Nellcor (1011) & Masimo(1049) with additional adapter. Call with MultiMed P/N and for additional information							
<b>Simed</b>								
S-100/ S-50/ S-100E Software Version E	1019	10'/3m	N/A	-	600333-00			
<b>Spacelabs</b>								
(Nellcor Based) 90351-D/G, 90465, 90303, 90466, 90467, 90489, 90651A-08, 90651A, F.M. IM77	1011	4'/1.2m	1020A	8'/2.4m	Sensor: (DS-100A) 690-0003-00 Adapter: 700-0002-00			
	1020N	10'/3m	N/A	-				
(Nova. Based Old Non S.B.) 90419, 90604, 90651-03	1020	10'/3m	1020B	1'/.3m	Sensor:(Nova Non-S.B.):369083-001 Adapter: 175-0646-00		 Female	
(Nellcor Based) All new modules (90651-08, 90496) Ultraview. 90651 updated modules with Ultraview Conn. Only	1011	4'/1.2m	1020E	8'/2.4m	Sensor: (DS-100A) 690-0003-00 Adapter: 700-0030-00			
	1020T	10'/3m	N/A	-				

	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Spacelabs</b>	<b>Continued</b>							
(Nellcor Based) Oximax technology	N/A	-	2003	8'/2.4m	700-0792-00			
(Nova. Based) All new modules (90651-08, 90496) Ultraview. Updated Newer Ultraview modules only w/ new 10-P Ultraview Conn.	1020S	10'/3m	1020D	1'/.3m	Sensor: (Nova S.B) 015- 0130-00 Adapter: 700-0029-00		 Female	
(Nova. Based) Opt. Extension Cable for all Novamatrix based monitoring (1020/1020S)	N/A	-	1020G	10'/3m	Old P/N: 4942-00 New P/N: 012-0073-00		 Female	
90309 (Nellcor Technology) PC SCOUT Monitor	1011	4'/1.2m	1020PC	8'/2.4m	Sensor (DS-100A):690- 0003-00			
90343-05 (Nellcor Technology)	1011	4'/1.2m	1020H	2'/.6m	Sensor (DS-100A):690- 0003-00 Adapter: 700-0014-00			
<b>Takoaka</b>								
All monitors w/ BCI Comp. sensors	1001	4'/1.2m	1056A	8'/2.4m				
	1056	10'/3m	N/A	-				
<b>VetOx / VetCo / Heska</b>								
All monitors w/ BCI Comp. sensors	<b>SEE BCI (1001)</b>							
<b>Vitalog</b>								
	<b>SEE CRITICARE (1002)</b>							

	BBS Code# Sensor	Sensor Length	BBS Code# Adapter/Ex tension Cable	Cable Length	OEM P/N's	Sensor Connector	Adapter Sensor Connector	Adapter Monitor Connector
<b>Weinmann</b>								
Oxycount Mini (BCI Compatible)	<b>SEE BCI (1001)</b>							
<b>Welch Allyn</b>								
(Nonin technology) (6)52SXX, (6)21SO/P, (6)22SO/P, (6)23SP	<b>SEE NONIN (1014)</b>							
(Nellcor technology) Protocol Propaq / Propaq LT, Atlas (6)52NXX, (6)22NO/P, (6)23NP	<b>SEE NELLCOR (1011)</b>							
VSM-300 Vital Signs Series Use Nellcor Pulse Oximetry. 53NXX REQUIRES OXIMAX	<b>SEE NELLCOR (1011)</b>  Newer Vital Signs Monitors use Nellcors OxiMax Sensors. The 1011 Sensor is not compatible with OxiMax Monitors. Call with Monitor Model number and information to determine if the monitor is OxiMax or 1011 compatible.							
Protocol Series. Protocol Series and Handheld Monitors (Nellcor AND Masimo Technology)	Protocol Series handhelds and Multiparameter Monitors use either Nellcor or Masimo Sensors. In addition, these monitors have 10 pin "ODU" adapter cables. Please call with Model #, specify Nellcor or Masimo and for more information.							
<b>ZOLL</b>								
M Series, AED Series	Masimo - N/A for repair only							